Solución GERMfree
Solution GERMfree
GERMfree solution

ADR: 3 F1 II UN 1219 IMDG: 3 II UN 1219 IATA: 3 II UN 1219 PAX: 353 CAO: 364
H225-H319-H336 P243-P261-P280-P304+P340-P305+P351+P338-P337+P313

SPECIFICATIONS:

APPLICATION:

Disinfection of surfaces, materials and equipment.

Antimicrobial properties against a wide spectrum of microorganisms: E.coli, Coliforms, Listeria, Staphylococcus aureus, Enterococcus faecalis, Bacillus subtilis, Candida albicans, MRSA, Pseudomonas aeruginosa, Streptococcus pyogenis, Proteus mirabilis, Presumptive clostridium perfringers, Presumptive enterobacteriaceae, fungi, salmonella, among others.

PHYSICO-CHEMICAL PROPERTIES:

Composition	Isopropanol 70%, water
Density (20°C)	8,49 g/cm ³
Storage conditions	Room temperature, avoid ignition sources nearby
Validity	24 months from manufacturing date (min. 6 months at delivery)
Appearance	Colourless liquid
Odour	Alcohol characteristic odour
Water solubility	Completely soluble
Self-ignition temperature	> 400 °C
Flammability	Highly flammable

DIRECTIONS OF USE:

Spray GERMfree directly on surfaces to disinfect until completely wet. Let the surface dry. It does not leave any residues.

PRESENTATION:

DSNF-GFS-750 GERMfree solution in spray, 750 ml

DSNF-GFR-750 GERMfree refill, 750 ml

Note:

Material Safety Data Sheets and Certificates of Analysis are downloadable free of charge at www.labbox.com. Disclaimer: labkem chemical products are exclusively for laboratory use and by instructed and experienced users. These chemical products may be dangerous, so please follow instructions mentioned on the Material Safety Data Sheets. Due to shipping restrictions or local/national regulations, all products may not be available in all countries or for some particular customers.

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SPECIFICATIONS:

APPLICATION:

Disinfection of surfaces, materials and equipment.

Anti-mycoplasma sp activity. It is also effective against Staphylococcus aureus, Streptococcus faecalis, Escherichia coli, Proteus mirabilis, Pseudomonas aeruginosa, Candida albicans, HBV / HIV, Papovavirus, Adenovirus, Poliovirus, among others.

PHYSICO-CHEMICAL PROPERTIES:

Composition	Isopropanol, Ethanol, Glutaraldehyde and Enrofloxacin
Density (20°C)	approx. 1 g/cm ³
Storage conditions	Room temperature, avoid ignition sources nearby
Validity	12 months from manufacturing date (min. 3 months at delivery)
Appearance	Colourless liquid
Odour	Alcohol characteristic odour
Water solubility	Soluble
Self-ignition temperature	approx. 26,5 °C (DIN51755)
pH	approx. 7

DIRECTIONS OF USE:

SURFACES: Spray MYCOPLASMAfree directly on surfaces to disinfect until completely wet. Wait 5 minutes and wipe the surfaces with a clean and dry paper towel. It does not leave any residues.

EQUIPMENT: Wet a paper towel with MYCOPLASMAfree and wipe the equipment. Remove after with a clean and dry paper towel.

WARNING: It may cause discoloration on sensitive surfaces. Test before use.

PRESENTATION:

DSNF-MYS-250 MYCOPLASMAfree solution in spray, 250 ml

DSNF-MYR-500 MYCOPLASMAfree refill, 500 ml

Note:

Material Safety Data Sheets and Certificates of Analysis are downloadable free of charge at www.labbox.com. Disclaimer: labkem chemical products are exclusively for laboratory use and by instructed and experienced users. These chemical products may be dangerous, so please follow instructions mentioned on the Material Safety Data Sheets. Due to shipping restrictions or local/national regulations, all products may not be available in all countries or for some particular customers.

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SPECIFICATIONS:

APPLICATION:

Decontamination/removal of DNA and RNA from surfaces, equipment, and pipettes from PCR workstations (POLYMERASE CHAIN REACTION).

PHYSICO-CHEMICAL PROPERTIES:

Composition	Phosphoric acid (< 5%), Anionic surfactants (< 5%)
Density (20°C)	approx. 1 g/cm ³
Storage conditions	Room temperature
Validity	12 months from manufacturing date (min. 3 months at delivery)
Appearance	Colourless liquid
Odour	Odourless
Water solubility	Soluble
Boiling point	> 100 °C
pH	approx. 1,3

DIRECTIONS OF USE:

ZeroDNA is a preparation designed to be used in PCR workstations. Particularly active against contamination by DNA and RNA in surfaces, laboratory equipment and pipettes.

Decontamination of surfaces (non porous): Spray ZeroDNA directly on the surface and spread with a piece of paper. Wash straight away with DNA-free water and dry with a clean paper towel.

Decontamination of laboratory devices and porous surfaces: Apply ZeroDNA in a paper towel and vigorously spread it in surfaces, devices or equipment. Wash straight away with DNA-free water and dry with a clean paper towel.

Decontaminate pipettes: Follow the manufacturer's instructions to disassemble pipettors. Dip pipette's rod, without seals, in ZeroDNA solution for one minute. Wash straight away with DNA-free water, and dry conveniently. Proceed to pipette reassembling.

It is recommended to spot test sensitive surfaces prior to use.

PRESENTATION:

DSNF-ZDS-250 ZeroDNA solution in spray, 250 ml

DSNF-ZDR-500 ZeroDNA refill, 500 ml

Note:

Material Safety Data Sheets and Certificates of Analysis are downloadable free of charge at www.labbox.com. Disclaimer: labkem chemical products are exclusively for laboratory use and by instructed and experienced users. These chemical products may be dangerous, so please follow instructions mentioned on the Material Safety Data Sheets. Due to shipping restrictions or local/national regulations, all products may not be available in all countries or for some particular customers.

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