**EVA180-B Rotary evaporator**

*Please read the User Manual carefully before use, and follow all operating and safety instructions!*

**user manual**

english / español

# User Manual

## LBX EVA180-B Rotary evaporator

### **Preface**

Users should read this manual carefully, follow the instructions and procedures, and beware of all the cautions when using this instrument.

### **Service**

If help is needed, you can always contact your dealer or Labbox via [www.labbox.com](http://www.labbox.com) (declare an incidence)

Please, provide the customer service representative with the following information:

- Serial number (on the back side)
- Description of the problem
- Your contact information

### **Warranty**

This instrument is guaranteed to be free from defects in materials and workmanship under normal use and service, for a period of 24 months from the date of invoice. The warranty is extended only to the original purchaser. It shall not apply to any product or parts which have been damaged on account of improper installation, improper connections, misuse, accident or abnormal conditions of operation.

For claim under the warranty, please contact your supplier.

## 1. Safety instructions

	Connect the equipment to a grounded power source to ensure the safety of the instrument and the experiment; connect power as per the equipment requirements.
	The use of this equipment in flammable, explosive, toxic or highly corrosive experiments is prohibited.
	Place the equipment on a horizontal, flat and stable table, creating a free space of 30 cm per side, allowing you to work comfortably. Make sure the surface is adequate for the weight of the equipment.
	The equipment must be used by qualified personnel who are familiar with the equipment and its operation through the user manual.
	Do not place the equipment in areas close to heat sources.
	While the machine is in operation, do not touch the heating surface to avoid burns.
	During operation, hazardous material such as flammable liquids or pathological material must be away from this area.
	Read the instruction manual before using the equipment and keep it on a place on which it can be accessed by anyone.
	Use the appropriate safety elements, such as eye protection, gloves, according to the information provided by the reagent's manufacturer.
	Pay special attention to splashes and protect yourself adequately.

-During the use of the equipment, use personal protective measures to avoid risks of possible damage such as:

- Splash burns and evaporation of liquids
- Poisoning by the emission of toxic or combustible gases
- Cuts or physical damage due to the use of fragile material such as glass.

-Beware of the dangers caused by:

- Flammable materials or media with a low boiling temperature
- Overfilling of the container
- Unsafe containers

- Make sure that the glass assembly has no stress due to the risk of breaking from causes as:
  - Mechanical hazards from the outside
  - Excess temperature in the equipment itself
- The temperature must be always set, at least, 25°C below the ignition temperature of the materials used
- Check equipment and accessories before each use and do not use components in poor condition
- Safe operation of the equipment is only guaranteed with the accessories described in chapter "Packing list". The accessories must be firmly attached to the equipment and cannot come loose
- Always turn off the power before installing accessories
- Do not cover the ventilation outlets of the unit to ensure adequate cooling of the rotary evaporator
- Adapt the quantity and type of solvent to the size of the distillation equipment
- Please know that high concentrations of gases, steam or particular matters can escape through the cooler. Avoid this risk by using, for example, cold water traps, gas wash bottles or an effective extraction system, such as showcases
- The glassware is designed for operations with maximum vacuum of 10 mbar. The glass must be cooled down before its removal and reuse
- Avoid formation of peroxide. Organic peroxides can accumulate in the distillation and generate exhaust residues that explode easily when decomposed
- Never heat the evaporation flask in the bath without first actuating the rotary motion. In this way, you will avoid possible losses of your sample
- Pay attention to the temperature setting; never leave the equipment unattended during operation and with the heating function on
- Do not handle the rotary evaporator when the rotating piston is working. It is recommended to gradually increase and/or decrease the rotation speed if situations such as an irregular operation occur
- The equipment can only be disconnected from the power outlet by pulling the plug base, not the cord
- The voltage indicated on the equipment must correspond to the electrical network used
- Make sure that the main power cord is not in contact with the heating surface. Do not cover the equipment
- Keep the equipment away from high magnetic fields

## 2. Rules

This equipment is specially designed for distillations for educational, laboratory or industrial purposes. This apparatus is not suitable for use at home or in environments that could be hazardous to the user or the equipment.

With the accessories recommended by the manufacturer, this equipment is capable of the following actions:

- Rapid distillation of liquids
- Evaporation of solutions and suspensions
- Crystallization, synthesis or cleaning of chemicals
- Recycle solvents

This instrument is intended for use in the following environments:

- Altitude:  $\pm 2,000$  m ASL
- Ambient temperature: 0-40°C
- Voltage fluctuation: within the range of -10%~+10% of normal value (the product is designed for indoor socket)

## 3. Inspection

### 3.1 Reception

Carefully unpack the instrument and check that the equipment and/or accessories have arrived undamaged. If necessary, contact the supplier who supplied the equipment for technical assistance.

**Note:**

- If the equipment is damaged, it must not be connected to an electrical source.
- Ensure that the device is unlocked before turning it on.

### 3.2 Packing list

The package includes the following items:

Item	Quantity
Main Unit	1
Heating bath	1
Vertical glass assembly	1
Power cable	1
User manual	1

Table 1

## 4. Control elements

### 4.1 Parts and controls


Fig. 1





Fig. 2

Fig. 3


Items	Descriptions
LED display	Indicates the current working status and parameter setting.
Forward indicator	The indicator is on when the instrument is in the state of forward rotation only
Reverse indicator	The indicator is on when the instrument is rotating forwardly and reversely in turn
Forward and reverse rotation interval setting button	Press this key to activate the setting of the forward and reverse rotation interval time. By turning the adjustment knob you can set the intermittent left and right rotation time within the range of 1-999 sec. When is set to 0 sec, the instrument rotates forward only
Adjustment knob	Turn the knob clockwise to increase the parameter and counter clockwise to decrease the parameter Press this key to start / stop the instrument
Water mode indicator	The indicator is on when the heating bath is in the water bath mode
Oil mode indicator	The indicator is on when the heating bath is in the oil bath mode
Heating indicator	The indicator is on when the Heating bath is in the heating state
Temperature adjustment knob	Turn the knob clockwise to increase the parameter and counter clockwise to decrease the parameter Press this key to start / stop the instrument Press the key for 3s to switch the heating mode
Heating bath power switch	Turn on / off the main power supply to the heating bath

**Table 2**

## 4.2 Display

After the main unit is turned on, the display is on and indicates the real time speed:

**Fig. 6**


This interface indicates the real-time speed. Turn the adjustment knob to set the speed and press the adjustment knob to start the instrument. The rotation speed can be changed by turning the adjustment knob while the instrument is in operation. After setting, the interface automatically changes to real-time speed. Press the “**INT**” key to stop the motor, with LED indicating “txxx” (xxx is the intermittent operation interval, to be set within 1-999sec), as shown in the figure.

**Fig. 7**

Turning the adjustment knob can adjust the forward and reverse rotation interval time. The instrument only rotates forward when the forward and reverse rotation interval time is 0. After setting, by pressing the “**INT**” key, the interface changes to the real-time speed display.

When the heating bath is turned on, water mode is the default mode, the water mode indicator is on and the interface is on, indicating the real-time temperature:


**Fig. 8**

Turn the temperature adjustment knob to set the target temperature of the heating bath. After setting, by pressing the knob, the heating indicator is on and the heating process begins. In the heating process, the target temperature can be changed by turning the knob. After setting, the interface automatically changes to real-time temperature.

Press the knob for 3s to switch the heating mode, with the corresponding indicator being on.


Water bath mode display:

**Fig. 9**

Oil bath mode display:

**Fig. 10**


## 5. Installation and operation

This chapter includes how to install the equipment and instructions for a first set-up.

**Note:** *Inspect the instrument before first use and during unpacking. If necessary, contact the supplier who supplied the equipment and keep the original packaging.*


### 5.1 Installation

Place the main unit on a flat table and squeeze the handle switch (a) to lift it to the highest position


**Fig. 11**

Tighten the screw (b) by hand, adjust the angle of the motor head and tighten the rotating head locking nut (c) to prevent accidental rotation.


**Fig.12**

Rotate the locking device (d) (Fig. 12) by 60° counter clockwise to mount the evaporating tube in place, and rotate the locking device (d) clockwise for locking

**Fig.13**

Mount the seal ring (Fig. 13 and Fig. 14)


**Fig.14**

Mount the condenser (f) (Fig. 15 and Fig. 16)

- Slide the locking nut (g) and pressure spring (h) onto the condenser.
- Screw down the locking nut (g) by hand


**Fig.15**

**Fig.16**

Mount the fixing device for the vertical condensation glass assembly (Fig. 17)

G

**Fig.17**

## Piping connection (Fig. 18)


Fig.18


Fig. 19

**Note:**

To ensure a good sealing performance, each glass ground joint may be coated with a layer of sealing grease.

## 5.2 Testing

- Check that voltage specified for this device matches with power network voltage
- Power socket should be grounded properly
- Turn on the power switch
- Turn the adjustment knob to set the target speed
- Press the forward and reverse rotation interval setting button to set the forward and reverse rotation interval time. If set to “0”, the instrument only rotates forward.
- Press the temperature adjustment knob for several seconds to switch between the water and oil bath modes.
- Turn the temperature adjustment knob to set the target temperature (up to 90°C under the water bath mode).
- Press the temperature adjustment knob to start heating
- Squeeze the handle switch to place the evaporating flask into the heating bath
- Press the adjustment knob to start rotation

**Note:**

Please don't remove the evaporating flask and receiving flask while the instrument is in operation.

If the vacuum pump is not used, keep the vacuum tube interface open.

## 6. Faults

This chapter helps to summarize the operations after a minor problem has occurred with the instrument. All possible problems are listed and also, their possible solutions.

### 6.1 General malfunction and its solutions

The equipment does not turn on:

-Check if the power line is disconnected

-Check that the cable and fuse are correctly connected to the equipment and to the electrical network

- ERR 001 - indicating driving motor error
  - Check the revolving shaft for any blockage and remove the foreign object if any.
- ERR 002 - indicating abnormal speed.
  - Turn off the instrument and restart it.
- ERR 003 - indicating excessively high temperature.
  - Turn off the heating bath and restart it when its temperature falls to the room temperature.
- ERR 004 - refers to a sensor failure
  - Change the temperature sensor

***If the above faults cannot be cleared, please contact your supplier.***

## 7. Maintenance and cleaning

Operate and maintain the product properly to keep a good working condition and extend the service life of the device. In routine service, keep the product dry and clean.

Remove the spilled liquid quickly on the external surfaces clean the outer surface with a mild cleaner, and do not connect the power supply until all surface is completely dry.

- If liquid or moist solid enters the product, please disconnect the power supply quickly and contact your supplier for more advice.
- If some liquid is spilled on the device surface, remove it quickly.

**WARNING*****Electrical hazard***

Perform the maintenance work with the equipment disconnected from the power supply and away from any flammable sources or steam.

Dye	Isopropanol
Building materials	Aqueous solution / isopropanol with active agent
Cosmetic	Aqueous solution / isopropanol with active agent
Food	Aqueous solution with active agent
Fuel oil	Aqueous solution with active agent

**Table 3**

## 8. Applicable regulations

<b>The product structure conforms to the following safety standard:</b>
EN 61010-1
UL 3101-1
CAN/CSA C22.2 (1010-1)
EN 61010-2-10
<b>The product structure conforms to the following EMC standard:</b>
EN 61326-1
<b>Conforming to the following EU standards:</b>
EMC Directive: 89/336/EWG
Instrument directive: 73/023/EWG

**Table 4**

## 9. Technical parameters

Items	Main unit	Heating bath
Voltage [VAC]	100-240	200-240/100-120
Frequency [Hz]	50/60	50/60
Power [W]	45	1010
Motor	DC brushless	-
Speed Range [rpm]	20-200	-
Speed Display	LED	LED
Temperature Range [°C]	RT -180	RT -180
Heat Control Accuracy [°C]	-	±1.5 (Oil ±3)
Lifting method	Manual	-
Lifting distance [mm]	110	-
Forward and reverse rotation interval setting	1-999	-
External Dimensions [D×W×H, mm]	440x320x450	300x300x240
Weight [kg]	7	3
Permissible Ambient Temperature [°C]	5-40	5-40
Permissible relative humidity	80%	80%
Protection Class acc. To DIN E N60529	IP20	IP20
USB	Yes	Yes

**Table 5**

## 10. Accessories and spare parts


This chapter includes accessories, sealing items and glassware that can be provided.

Use only the accessories of LABBOX to guarantee the operation of the equipment and ensure a good functioning of the system.

### 10.1 Glass assembly

Labbox reference	Description
REVA-A01-001	Vertical glass set, includes condenser, distillation flask (269/32) and collector flask (35/20)
REVA-A02-001	Collecting flask (35/20), capacity 1000ml
REVA-A03-001	Piece steam step
REVA-A04-001	Vacuum seal

**Table 6**

**Fig. 20**

\*


Evaporation flask	Capacity
FLEG-051-001	50 ml
FLEG-101-001	100 ml
FLEG-250-001	250 ml
FLEG-500-001	500 ml
FLEG-1K0-001	1000 ml
FLEG-2K0-001	2000 ml

**Table 7**


## 10.2 Vacuum references and accessories

Labbox reference	Description
REVA-A04-001	Vacuum seal
REVA-A05-001	Evaporating flask clamp
LTXT-007-005	Latex tubing, amber, for 8/9 hose, 5m
RUTT-008-002	Rubber vacuum tubing for 9/10 hose, 2m
VACU-R20-001	Vacuum pump LBX R20 chemically resistant, 18 l/min, -750 mm Hg

**Table 8**

**Fig. 21**

For more information about accessories or laboratory equipment that you need, visit our website [www.labbox.com](http://www.labbox.com) or call +0034 902 99 53 18.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

## 11. Storage and transportation

This chapter introduces you on how to pack away the instrument, how to prepare it for transport, specific storage and transport conditions.

	<p><b>Alert</b> <b>Biohazard</b></p> <ul style="list-style-type: none"><li>- If you require maintenance service, the instrument must be cleaned before shipment to avoid pollution of hazardous substances and also to be sent back in the original packaging</li><li>-If the instrument will not be used for a long time, please, switch off and store it in a dry and clean place with a steady room temperature.</li><li>-Ensure that the device will not get wet or hit during transport</li></ul>
---	--

	<p><b>Alert</b> <i>Electrical hazard</i></p> <p>Always pull out the socket to avoid having electric current nearby</p>
---	--

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