
EASY

Bottle top dispenser EASY 20K+

Please read the User Manual carefully before use, and follow all operating and safety instructions!

user manual

english / español

User manual

Bottle top dispenser, EASY 20K+

Introduction

Users must read this manual carefully before using the equipment and follow the instructions and procedures in order to be aware of all the precautions, as well as in order to obtain the maximum benefits and duration of the equipment.

Service

If help is needed, contact with your supplier/manufacturer or Labbox via:

www.labbox.com

Please, provide the customer service representative with the following information:

- Serial number (on the back side)
- Description of the problem
- Your contact information

Warranty


This instrument is guaranteed to be free from defects in materials and workmanship under normal use and service, for a period of 24 months from the date of invoice. The warranty is extended only to the original purchaser. It shall not apply to any product or parts which have been damaged on account of an improper installation, improper connections, misuse, accident, or abnormal conditions of operation. For claim under the warranty, please contact your supplier.

1. Safety instructions

	Place the Bottle top dispenser EASY 20K on a flat, horizontal, and stable table, leaving at least 30 cm of free space on each side. Make sure that the surface is adequate for the device's weight.
	This item must be used only by previously qualified staff that is familiar with the equipment and knows how to operate it according to the instructions manual.
	Carefully read the instructions in this manual before using the device.
	Use the necessary protective clothes, eye protection and gloves according to the reagent's manufacturer's information.
	Pay special attention to splashes and protect yourself adequately

- Do not operate this instrument in any manner not described in this 'User manual'.
- In case of abnormal circumstances during the manipulation of this device, please stop the operation of the equipment and contact your supplier.
- The precautions described in this 'User manual' are carefully developed in an attempt to cover all the possible risks. However, it is also important that you are alert for unexpected incidents.
- Use this instrument only for dispensing liquids, with strict regard to the defined operating exclusions and limitations. (chart on p.11)
- Explosive liquids are also forbidden to use with this dispenser.
- Before use, check the instrument for visual damages (piston difficult to move, sticking valves or leakage...) in case of trouble, immediately stop dispensing.

2. Components

4.1 Bottle top dispenser

1. Piston cover
2. Piston seat
3. Cylinder sleeve
4. Volume adjustment knob
5. Piston (inside the socket)
6. Holding grip
7. Air vent cap
8. Bottle top screw seat
9. Filling valve assembly
10. Filling tube
11. Discharge tube
12. Discharge tube sleeve
13. Recovery switch
14. Closure cap
15. Recovery valve
16. Seal

3. Application and design

This instrument is designed for dispensing liquids, observing the following limits:

- 15 to 40 °C of instrument and reagent
- Vapor pressure up to 1000 mbar
- Density up to 2.2 g/cm³

4. Materials

When the instrument is correctly used, the dispensed liquid comes into contact with only the following chemically resistant materials: Borosilicate glass, PFA, FEP, PTFE, PVDF, ETFE, Hastelloy; PP (closure cap).

5. Operating exclusions

Never use this instrument for

- liquids attacking FEP, PFA, PTFE or Hastelloy.
- liquids attacking borosilicate glass (e.g., hydrofluoric acid)
- liquids which are decomposed catalytically by Hastelloy (e.g., H₂O₂)
- explosive liquids
- suspensions of solid particles which may clog or damage the instrument.

6. Operating limitations

Concentrated hydrochloric acid and concentrated nitric acid, chlorinated and fluorinated hydrocarbons, or liquids which form deposits may make the piston difficult to move or may cause jamming.

Attention:


If there is a sign of a potential malfunction (e.g., piston difficult to move) never use force. Immediately stop dispensing and follow cleaning instructions or contact the manufacturer.

Note:

Check the compatibility table or contact the manufacturer for special applications (e. g., trace material analysis).

7. Use

7.1. Connecting the filling tube


Push the filling tube into the valve block as deep as possible.


7.2. Mounting the instrument on the bottle

The instrument can be screwed directly onto a GL 45 screw thread. For other bottles, use the accompanying adapters.

7.3. Priming the instrument

- Hold the discharge tube and carefully remove the closure cap.
- Move the closure cap backwards along the discharge tube, away from the discharge tube opening.


- Hold an appropriate collecting vessel under the discharge tube opening.
- Gently pull the piston up, then press it down rapidly.
- Repeat this procedure until most of the air bubbles in the glass cylinder and discharge tube have been removed.

Note:

Before using the instrument for the first time, ensure it is rinsed carefully or discard first few samples dispensed.


8. Priming/Dispensing

Attention:

Follow all Safety Instructions and observe the Operating Exclusions and Limitations. The discharge tube must always point away from the user or any other person. Never press down the piston when the closure cap is pushed on. Reagent may drop off the discharge tube and closure cap.

8.1. Setting the volume

- Rotate the volume adjustment knob counterclockwise to unlock it
- Move the knob along the scale to target volume
- Rotate the knob clockwise to lock it at the target volume.


8.2. Fill and dispense

- Hold an appropriate collecting vessel under the discharge tube opening.
- Pull the piston up slowly and evenly to the upper stop to draw reagent into the glass cylinder.
- Push the piston down gently and evenly to dispense the reagent into the vessel.


9. Cleaning

The instrument must be cleaned in the following situations to assure correct operation:

- if the piston becomes difficult to move
- after dispensing liquids which form deposits
- before changing the reagent
- prior to long term storage
- should liquid have accumulated in the closure cap
- prior to sterilization
- before replacing the valves

1. Remember to recover the reagent before cleaning the dispenser.
2. Make sure the dispenser is still properly mounted on the reagent bottle and there is plenty space in the bottle for the reagent inside the dispenser to be recovered.
3. Rotate the recovery switch clockwise by 90 degrees to the direction perpendicular to the discharge tube, and the remaining reagent inside the dispenser will flow back into the bottle underneath.
4. Turn off the recovery switch. Carefully remove the dispenser from the reagent bottle to another empty bottle. Pull and push the piston for several time to dispense the remaining reagent from the discharge tube into the reagent bottle as much as possible.
5. Remove the dispenser from the empty bottle to another bottle full of distilled water (or other proper wash buffer).
6. Rinse the instrument by dispensing the washing liquid into a beak again and again.
7. Empty the bottle, flush it and fill it with distilled water. Repeat step 6 to 7.

8. Empty the instrument.

9. Unscrew the piston seat (1).

10. Carefully pull the piston out of the glass cylinder by pulling on the piston seat only.

11. Carefully remove sediments at the edge of the glass cylinder with a screwdriver. Clean piston and cylinder with a bottle brush.

12. Carefully reassemble the instrument in the reverse order. Flush it with distilled water.

Note:


Insert the piston straight into the cylinder, not at an angle.

9.1. Clean/Replace the valve

1. Follow all cleaning instructions.

2. Remove the filling tube.

3. Unscrew the filling valve using the mounting tool. Remove the valve with its seal. Make sure the seal does not remain in the valve block.

4. Rinse the filling valve in the cleaning solution, and clean with soft brush. Should valve ball become sticky, push it into the filling opening of valve with a pointed instrument, and release it.
5. Screw the cleaned/new filling valve together with its seal into the valve block and fasten it firmly using the mounting tool.

10. Autoclaving

This instrument withstands steam sterilization at 121 °C, 1 bar absolute (15 psi) for 20 minutes in accordance with general guidelines. It is the user's responsibility to ensure effective autoclaving.

10.1. Preparation prior to the autoclaving

1. Clean the instrument before sterilization.
2. Remove the filling tube (10) and put it on a towel.
Avoid contact with hot metal surfaces.
3. Lay the dispenser with piston at the bottom extreme position on the towel.
4. Autoclave the dispenser.

Note:

Wait the dispenser to cool down in the closed autoclave to avoid deformation through too rapid cooling. After every autoclaving, inspect all parts for deformities or damage. If necessary, replace them. Do not reassemble the instrument until it has cooled down to room temperature (cooling time approx. 2 hours).

Mounting tool can be steam sterilized at 121 °C.

11. Checking the volume

The accuracy and coefficient of variation of the instrument are determined gravimetrically as follows:

- Set to the nominal volume.
- Dispense distilled H₂O.
- Weigh the dispensed quantity on an analytical balance.
- Calculate the dispensed volume taking the temperature into account.
- Perform at least 10 dispensing and weighing operations.
- Calculate the accuracy (A%) and the coefficient of variation (CV%) by means of the formular used in statistical quality control. The proceeding is described e.g., in DIN EN ISO 8655-6. Observe the operating manual of the balance manufacturer and the corresponding standards.

Mean value: $\bar{x} = \Sigma x_i / n$,

$x_i = results\ of\ weighings$

$n = number\ of\ weighings$

Mean volume: $\bar{V} = \bar{x} \cdot z$

$z = correction\ factor\ (e.g.\ 1.0029\ \mu l/mg\ at\ 20^\circ C, 1013hPa)$

Accuracy: $A\% = (\bar{V} - V_0 / V_0) \cdot 100$,

$V_0 = nominal\ volume$

Standard deviation: $s = Z \cdot \sqrt{\Sigma (x_i - \bar{x})^2 / n - 1}$

Coefficient of variation: $CV\% = 100 / \bar{V}$

12. Trouble shooting

Problem	Possible cause	Corrective action
Piston with difficulty to move	Dirt, possible formation of crystals	Stop dispensing immediately. Loosen the piston by circular motion without disassembling. Follow all the cleaning instructions.
Filling not possible	-Volume adjusted to minimum setting -Sticking filling valve	-Set the required volume -Clean the filling valve and replace with sealing washer if necessary.
Instrument does not dispense liquid	- Closing cap is on -Discharge tube not firmly connected or damaged	-Remove closing cap -Push on the closing cap. Hold the piston in place and push the original manufacturer's discharge tube on firmly. Replace a deformed or damaged discharge tube.
Air bubbles in the instrument	- Reagent with high vapor pressure has been poured too quickly. -Seal not inserted, knurled locking nuts not firmly connected -Air not expelled from the instrument -Filling tube not firmly connected or damaged -Valves dirty, not firmly installed, or damaged	-Pour the reagent slowly. -Double check that seal (12) is in place, and that the knurled locking nuts at the discharge tube are firmly seated and finger tightened -Priming the instrument -Push the filling tube on firmly. If necessary, cut off approx. 1 cm of tube at the upper end and reconnect it or replace filling tube. -Cleaning procedure. Tighten the valves using the mounting tool. If necessary, replace the valves and sealing washers.
Dispensed volume is too low	-Discharge tube not firmly connected or damaged -Filling tube not firmly connected or damaged -Valves not firmly installed or damaged -Sticky filling valve	-Push the original manufacturer's tube on firmly. -Cleaning procedure. Push the filling tube on firmly. If this does not suffice, cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube. -Clean the instrument. Tighten the valves using the mounting tool, if necessary, replace the valves and sealing washers. - Clean or replace discharge tube, if necessary.

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