

GENERAL INFORMATION

Product name : Spatula, Chattaway**Description :** Made of stainless steel

TECHNICAL DATA

reference	l(mm)	wspatula (mm)	pcs/pack
SPCS-150-005	150	6	5
SPCS-200-005	200	7	5

PACKAGING AND LOGISTICS

reference	vol (l)	kg	TARIC	GTIN
SPCS-150-005	0,0014	0,055	82055980	08434868043947
SPCS-200-005	0,08575	0,075	82055980	08434868043954

PRODUCT PHOTO



INFORMACIÓN GENERAL

Nombre del producto : Espátula Chattaway**Descripción :** De acero inoxidable

DATOS TÉCNICOS

referencia	l (mm)	wespátula (mm)	unidades por ref.
SPCS-150-005	150	6	5
SPCS-200-005	200	7	5

EMBALAJE Y DATOS LOGÍSTICOS

referencia	vol (l)	kg	TARIC	GTIN
SPCS-150-005	0,0014	0,055	82055980	08434868043947
SPCS-200-005	0,08575	0,075	82055980	08434868043954

FOTO DEL PRODUCTO



INFORMATIONS GÉNÉRALES

Nom produit : Spatule Chattaway**Description :** En acier inoxydable

DONNÉES TECHNIQUES

référence	long. (mm)	larg.spatule (mm)	unités par ref.
SPCS-150-005	150	6	5
SPCS-200-005	200	7	5

EMBALLAGE ET LOGISTIQUE

référence	vol (l)	kg	TARIC	GTIN
SPCS-150-005	0,0014	0,055	82055980	08434868043947
SPCS-200-005	0,08575	0,075	82055980	08434868043954

PHOTO PRODUIT



INFORMAZIONE GENERALE

Nome del prodotto : Spatola Chattaway**Descrizione :** In acciaio inossidabile

DATI TECNICI

referenza	L (mm)	hspatola (mm)	unità per ref.
SPCS-150-005	150	6	5
SPCS-200-005	200	7	5

IMBALLAGGIO E DATI LOGISTICI

referenza	vol (l)	kg	TARIC	GTIN
SPCS-150-005	0,0014	0,055	82055980	08434868043947
SPCS-200-005	0,08575	0,075	82055980	08434868043954

FOTO DEL PRODOTTO

