

Tabla de compatibilidad

Embudo Büchner + junta de goma Guko + matraz Kitasato

	Referencia	EMBUDO BÜCHNER DE PORCELANA PREMIUM LINE									EMBUDO BÜCHNER DE PP		EMBUDOS SEGÚN DR. HIRSCH PREMIUM LINE			
		BFUK-035-001	BFUK-100-001	BFUK-200-001	BFUK-350-001	BFUK-600-001	BFUK-800-001	BFUK-1K1-001	BFUK-2K0-001	BFUK-4K8-001	BFUP-070-001	BFUP-110-001	HFFK-030-001	HFFK-040-001	HFFK-050-001	HFFK-075-001
MATRAZ KITASATO, FORMA ERLLENMEYER (GLASSCO)	FFE3-100-002	GUKO-002-001	GUKO-003-001	GUKO-003-001	-	-	-	-	-	-	GUKO-003-001	-	GUKO-002-001	GUKO-002-001 GUKO-003-001	GUKO-002-001 GUKO-003-001	-
	FFE3-125-002	GUKO-003-001	GUKO-003-001	GUKO-003-001	-	-	-	-	-	-	GUKO-003-001	-	(GUKO-003-001)	GUKO-003-001	GUKO-003-001	-
	FFE3-250-002	GUKO-003-001 GUKO-004-001	GUKO-003-001 GUKO-004-001	GUKO-003-001 GUKO-004-001	GUKO-004-001	-	-	-	-	-	GUKO-003-001 GUKO-004-001	-	-	-	GUKO-003-001 GUKO-004-001	GUKO-003-001 GUKO-004-001
	FFE3-500-002	GUKO-004-001	GUKO-004-001	GUKO-004-001	GUKO-004-001	GUKO-004-001	-	-	-	-	GUKO-004-001	-	-	-	(GUKO-004-001)	GUKO-004-001
	FFE3-1K0-002	GUKO-004-001	GUKO-004-001	GUKO-004-001 GUKO-005-001	GUKO-004-001 GUKO-005-001	GUKO-004-001 GUKO-005-001	GUKO-005-001	-	-	-	GUKO-004-001 GUKO-005-001	-	-	-	(GUKO-004-001)	GUKO-004-001 GUKO-005-001
	FFE3-2K0-001	-	GUKO-005-001	GUKO-005-001	GUKO-005-001	GUKO-005-001	GUKO-005-001	GUKO-005-001	-	-	GUKO-005-001	GUKO-005-001	-	-	-	(GUKO-005-001)
	FFE3-5K0-002	-	-	GUKO-006-001	GUKO-006-001	GUKO-006-001	GUKO-006-001	GUKO-006-001	GUKO-006-001	GUKO-006-001	-	GUKO-006-001	-	-	-	(GUKO-006-001)
MATRAZ KITASATO, FORMA ERLLENMEYER (LBG)	FFEN-250-001	GUKO-004-001	GUKO-004-001	GUKO-004-001	GUKO-004-001	-	-	-	-	-	GUKO-003-001 GUKO-004-001	-	-	-	GUKO-003-001 GUKO-004-001	GUKO-003-001 GUKO-004-001
	FFEN-500-001	GUKO-004-001	GUKO-004-001	GUKO-004-001	GUKO-004-001	GUKO-004-001	-	-	-	-	GUKO-004-001	-	-	-	(GUKO-004-001)	GUKO-004-001
	FFEN-1K0-001	GUKO-004-001	GUKO-004-001	GUKO-004-001	GUKO-004-001	GUKO-004-001	-	-	-	-	GUKO-004-001	-	-	-	-	GUKO-004-001
MATRAZ KITASATO, FORMA ERLLENMEYER DE PP	FFEP-250-001	GUKO-003-001	GUKO-003-001	GUKO-003-001 GUKO-004-001	(GUKO-004-001)	-	-	-	-	-	GUKO-003-001 GUKO-004-001	-	-	GUKO-003-001	GUKO-003-001 GUKO-004-001	GUKO-003-001 GUKO-004-001
	FFEP-500-001	GUKO-004-001	GUKO-004-001	GUKO-004-001	GUKO-004-001	GUKO-004-001	-	-	-	-	GUKO-004-001	-	-	(GUKO-004-001)	GUKO-004-001	GUKO-004-001
	FFEP-1K0-001	-	(GUKO-005-001)	GUKO-005-001	GUKO-005-001	GUKO-005-001	GUKO-005-001	(GUKO-005-001)	-	-	GUKO-005-001	GUKO-005-001	-	-	-	GUKO-005-001
	FFEP-2K0-001	-	-	GUKO-005-001	GUKO-005-001	GUKO-005-001	GUKO-005-001	GUKO-005-001	-	-	GUKO-005-001	GUKO-005-001	-	-	-	GUKO-005-001